

## Variable guide: Hepatitis B viral load

Please refer to the below table when completing ISOSS hepatitis B viral load test results. ISOSS requires hepatitis B viral load in IU/ml. If viral load is given in another format, please use the table to convert into IU/ml.

This is the format you may receive results in from your local laboratories

**Preferred format for ISOSS Online**

IU/ml	IU/ml in scientific notation	Log <sub>10</sub> – find IU/ml by typing '=10^X' in excel.
100 IU/ml	1.0 X 10 <sup>2</sup>	2.0
1,000 IU/ml	1.0 X 10 <sup>3</sup>	3.0
10,000 IU/ml	1.0 X 10 <sup>4</sup>	4.0
53,000 IU/ml	5.3 X 10 <sup>4</sup>	4.72
100,000 IU/ml	1.0 X 10 <sup>5</sup>	5.0
1,000,000 IU/ml*	1.0 X 10 <sup>6</sup> *	6.0*
2,000,000 IU/ml	2.0 X 10 <sup>6</sup>	6.301

\*Above these values are higher infectivity risk